



7NG 35GHLP 02 905



7NG 35NHPP 02 905



7NG 35GBLR 00 905



7NG 35GDLR 02 512




7NG 35NBPU 00 512

## POLYURETHANE ADJUSTABLE CHAIR

### PRODUCT ADVANTAGES :

- Articulated nylon pads or casters
- Height adjustment seat by chromed or black gas lift column, with shock absorber
- 4 available heights, for adjustments between 410 to 970mm. 130 and 240mm strokes
- Articulated backrest in height and depth
- Available polyamid footing with chromed branches, adjustable in height
- Black or blue flexible polyurethane foam seat and backrest, offering optimal comfort (4cm thickness)
- Does not overheat or cause perspiration  
Easy to clean using common household products

## CREATE YOUR KANGO CHAIR :



1	MATERIAL	4	MECHANISM	7	HEIGHT	10	SEATING OPTIONS	12	COLOR
7	Polyurethane	3	Adjustable chair	T	Very high ( <i>chromed</i> )	0	Without	512	Blue « roy »
				H	High ( <i>chromed or black</i> )			905	Black
2	FINISH	5	VERSION	D	Intermediate ( <i>chromed</i> )	11	FOOTREST		
N	Normal/Standard	5	Variant 3	B	Low ( <i>chromed or black</i> )	0	Without		
3	TYPE OF SEATING	6	COLUMN	8	BASE	2	Polyamide with chromed ring		
G	New PU range	G	Chromed gas lift column	L	Chromed long				
		N	Black gas lift column	P	Polyamid long				
9	GROUND CONNECTIONS								
P	Standard pads								
R	Standard casters								
V	Braked vacuum casters								
C	Locked under load casters								
U	Casters with soft rubber tread								



5NN 35GHLP 02 000



5NN 35NHPP 00 000



5NN 35GBLR 00 000



5NN 35NBPR 00 000

5NN 35NHPP 02 000

## WOODEN ADJUSTABLE CHAIR

### PRODUCT ADVANTAGES :

- Chromed 5 branch base, 590mm diameter or polyamid 5 branch base, 600mm diameter
- Articulated nylon pads or casters
- Height adjustment seat by chromed or black gas lift column, with shock absorber
- 4 available heights, for adjustments between 370 to 930mm. 130 and 240mm strokes
- Available polyamid footing with chromed branches, adjustable in height
- Multi-ply varnished beech wood seat and backrest, 10mm thickness
- Articulated backrest in height and depth
- Ecological and easy to clean

### CREATE YOUR KANGO CHAIR :

5
N
N
3
5
?
?
?
?
0
?
0
0
0

1 2 3

4 5 6 7 8 9

10 11

12

<b>1</b>	<b>MATERIAL</b>	<b>4</b>	<b>MECHANISM</b>	<b>7</b>	<b>HEIGHT</b>	<b>10</b>	<b>SEATING OPTIONS</b>	<b>12</b>	<b>COLOR</b>
5	Wood	3	Adjustable chair	T	Very high ( <i>chromed</i> )	0	Without	000	No colour
				H	High ( <i>chromed or black</i> )				
<b>2</b>	<b>FINISH</b>	<b>5</b>	<b>VERSION</b>	<b>8</b>	<b>BASE</b>	<b>11</b>	<b>FOOTREST</b>		
N	Normal/Standard	5	Variant 3	D	Intermediate ( <i>chromed</i> )	0	Without		
				B	Low ( <i>chromed or black</i> )	2	Polyamide with chromed ring		
<b>3</b>	<b>TYPE OF SEATING</b>	<b>6</b>	<b>COLUMN</b>	<b>9</b>	<b>GROUND CONNECTIONS</b>				
N	Normal/Standard	G	Chromed gas lift column	P	Standard pads				
		N	Black gas lift column	R	Standard casters				
				V	Braked vacuum casters				
				C	Locked under load casters				
				U	Casters with soft rubber tread				



6NN 35GBLR 00 705



4NN 35NBPR 00 905

6NN 35NHPP 02 705




6NN 35GHLP 02 505

## FABRIC/VINYL ADJUSTABLE CHAIR

### PRODUCT ADVANTAGES :

- Chromed 5 branch base, 590mm diameter or polyamid 5 branch base, 600mm diameter
- Articulated nylon pads or casters
- Height adjustment seat by chromed or black gas lift column, with shock absorber
- 4 available heights, for adjustments between 400 to 960mm. 130 and 240mm strokes
- Available polyamid footing with chromed branches, adjustable in height
- Seat and backrest, available in blue or grey fabric, or black vinyl, on high density foam with protective casing
- Articulated backrest adjustable in height and depth

### CREATE YOUR KANGO CHAIR :



<b>1</b>	<b>MATERIAL</b>	<b>4</b>	<b>MECHANISM</b>	<b>7</b>	<b>HEIGHT</b>	<b>10</b>	<b>SEATING OPTIONS</b>	<b>12</b>	<b>COLOR</b>
	4 Vinyl 6 Fabric		3 Adjustable chair		T Very high ( <i>chromed</i> ) H High ( <i>chromed or black</i> ) D Intermediate ( <i>chromed</i> ) B Low ( <i>chromed or black</i> )		0 Without		505 Dark blue 705 Dark grey 905 Black
<b>2</b>	<b>FINISH</b>	<b>5</b>	<b>VERSION</b>	<b>8</b>	<b>BASE</b>	<b>11</b>	<b>FOOTREST</b>		
	N Normal/Standard		5 Variant 3		L Chromed long P Polyamid long		0 Without 2 Polyamide with chromed ring		
<b>3</b>	<b>TYPE OF SEATING</b>	<b>6</b>	<b>COLUMN</b>	<b>9</b>	<b>GROUND CONNECTIONS</b>				
	N Normal/Standard		G Chromed gas lift column N Black gas lift column		P Standard pads R Standard casters V Braked vacuum casters C Locked under load casters U Casters with soft rubber tread				



4DL 31NHPR 02 705



4DL 31NBPP 00 505



6XL 31NHPP 02 505




6XL 31NBPR 00 905

## WIDE FABRIC/VINYL ADJUSTABLE CHAIR

### PRODUCT ADVANTAGES :

- Polyamid 5 branch base, 600mm diameter
- Articulated nylon pads or casters
- Height adjustment seat by black gas lift column, with shock absorber
- 2 available heights, for adjustments between 430 to 810mm. 130 and 240mm strokes
- Available polyamid footring with chromed branches, adjustable in height
- Wide seat and high backrest, available in black or blue M1 fabric, or in grey or blue M1 vinyl, on high density polyether foam
- Articulated backrest adjustable in height and depth

## CREATE YOUR KANGO CHAIR :



<b>1</b>	<b>MATERIAL</b>	<b>4</b>	<b>MECHANISM</b>	<b>7</b>	<b>HEIGHT</b>	<b>10</b>	<b>SEATING OPTIONS</b>	<b>12</b>	<b>COLOR</b>
4	Vinyl	3	Adjustable chair	H	High	0	Without	505	Dark/Cobalt blue
6	Fabric			B	Low			705	Dark/Elephant grey
<b>2</b>	<b>FINISH</b>	<b>5</b>	<b>VERSION</b>	<b>8</b>	<b>BASE</b>	<b>11</b>	<b>FOOTREST</b>	905	Black
X	Xtreme or Madrid	1	Variant 1	P	Polyamide long	0	Without		
D	Diabolo club	<b>6</b>	<b>COLUMN</b>			2	Polyamide with chromed ring		
G	Ginkgo	N	Black gas lift column	<b>9</b>	<b>GROUND CONNECTIONS</b>				
<b>3</b>	<b>TYPE OF SEATING</b>			P	Standard pads				
L	Wide			R	Standard casters				
				U	Casters with soft rubber tread				