

SciQuip Multi-purpose Syringe Diluter Dispenser

- 25ul – 25ml range
- 0.2% accuracy (0.002% of the nominal syringe volume)
- 4 operating modes when used with 2 syringes
- High torque precision stepper motor offers unrivalled positional accuracy
- 9 stage adjustable pipetting speeds
- 2 year warranty

Part number	Description
SP3060	SciQuip Multi-purpose 25uL-25mL Syringe Diluter Dispenser - Complete
SP3065	Syringe - 250ul
SP3066	Syringe - 500ul
SP3067	Syringe - 1ml
SP3068	Syringe - 2.5ml
SP3069	Syringe - 5ml
SP3070	Syringe - 10ml
SP3071	Syringe - 25ml



Technical Specifications

Accuracy	±1%
Precision	±0.2%
Flow Rate	0.003-6000 uL/s (depending on syringe)
Compatible Syringes	25, 50, 100, 250, 500uL, 1, 2.5, 5, 10, 25ml
Fluid Path Material	Borosilicate, PTFE, CTFE
Power requirements	100-240V, 50/60Hz,
Dimensions (WDH)	190x180x350mm
Weight	4.5kg

Multi-purpose

Dual syringe capability means that this unit can be used for dedicated dispensing / diluting as well as a dual syringe diluter, single syringe dispenser, dual syringe dispenser and a continuous dispenser.



SciQuip Ltd
Newtown
Wem
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk