

***Vortex VariMix & FixMix***





## SciQuip Vortex VariMix & FixMix

- Vortex shakers with compact design
- Touch or continuous operation
- Wide speed range 0 - 2,500rpm (VariMix)
- Multiple applications with additional accessories, microtubes, test tubes, microplates
- Robust case with specially designed silicone feet to absorb vibration at high speeds
- Heavy base with ultra high vibration damping

### Technical Data

Shaking movement	Orbital
Orbital diameter (mm)	4
Speed range (rpm)	0 - 2,500 (VariMix)
Operation	Touch or Continuous
Dimensions (w x d x h, mm)	127 x 130 x 160
Weight (kg)	3.5
Voltage (V)	200 - 240
Frequency (Hz)	50/60
Permissible ambient temperature (°C)	5 - 40
Permissible relative humidity (%)	80

### Vortex Mixers

Part No.	Description	Speed Range
SP2260-VM	Vortex VariMix (includes standard adapter)	Variable
SP2260-FM	Vortex FixMix (includes standard adapter)	Fixed

### Accessories

Part No.	Description	Operation
SP2260-VT1.1	Standard touch pad for tubes up to 30mm diameter	Touch/ Continuous
SP2260-VT1.3	Universal attachment plate	Continuous

### Accessories for universal attachment plate SP2260-VT1.3

Part No.	Description	Operation
SP2260-VT3.1	For 48 x microtubes	Continuous
SP2260-VT3.2	For 15 x 10mm diameter tubes	Continuous
SP2260-VT3.3	For 12 x 12mm diameter tubes	Continuous
SP2260-VT3.4	For 8 x 16mm diameter tubes	Continuous
SP2260-VT3.5	For 8 x 20mm diameter tubes	Continuous
SP2260-VT3.6	For 1 x Erlenmeyer flask 100 - 250ml	Continuous
SP2260-VT3.7	For 1 x microplate	Continuous



SciQuip Ltd  
Newtown  
Wem  
Shrewsbury  
Shropshire  
SY4 5NU, UK

Tel: +44 (0)1939 234222  
Fax: +44 (0)1939 234221  
Email: [info@sciquip.co.uk](mailto:info@sciquip.co.uk)  
[www.sciquip.co.uk](http://www.sciquip.co.uk)