

SciSpin

Centrifuges

SciSpin Microfuges



www.sciquip.co.uk

SciSpin Micro

- Compact design
- Maximum 15,000rpm/15,100g
- Capacity 12 x 1.5/2.2ml
- User friendly large LCD display
- Unique design minimises chamber temp
- Brushless motor
- Over speed & over temp detector
- Run-end alarm, pulse button, automatic lid open

Ordering Information

SciSpin Microfuge

Part No.	Description
SQ-6000	SciSpin Micro including 12-place rotor with snap on lid

Accessories

Part No.	Description
SQ-6010	Rotor for 12 x 1.5/2.0ml with snap on lid
SQ-6020	0.5ml adapter for use with rotor
SQ-6025	0.2ml adapter for use with rotor

SciSpin Micro

Capacity	12 x 1.5/2.2ml
Max Speed	15,000rpm (increment 100rpm)
Min speed	500rpm
Max RCF	15,100g (increment 100g)
Timer	30 sec, 99 min or continuous
Dimensions (WxDxH)	255 x 245 x 140mm
Weight	6kg
Noise level	56dB
Power	220V 50/60Hz



SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk